

Improved physical function using a power enhancing glove in persons with Inclusion Body Myositis

Objective: To investigate if a power enhancing glove is feasible to use in persons with IBM.

Method: Data were collected during The Myositis Association’s (TMA) patient conference in 2023. The study included 40 individuals with IBM self-reported hand weakness.

The participants rated their ability to perform self-selected tasks that were difficult for them to perform with and without Carbonhand. Participants selected three activities experienced as difficult and important from the modified Patient Reported Outcome-Upper Extremity Function Scale for IBM (PRO-UEFS). They performed the selected activities rating their perceived limitation on a 5-point scale (0 = “unable to do”, 4= “without any difficulty”) first without the Carbonhand glove and then upon using the glove (see table 1).

Results: The three most selected activities were lifting free weights, open previously opened jars, lift a heavy bag from the floor, holding av frying pan and picking a coin from a table. In the open-ended questions, participants documented that the glove would be beneficial for use in everyday tasks, lifting objects, grocery shopping, stabilizing the hand and would increase independence.

All activities were perceived easier to perform with the glove with ratings on perceived limitation reflecting “with little difficulty” (3) to “without difficulty” (4).

Chosen activity	Without the glove	With the glove
Lift a heavy bag from the floor n=15	2 (2-2)	4 (4-4)
Open previously open jars, n=16	1 (1-2)	4 (3-4)
Lifting free weights, n=16	2 (1-2)	4 (3.25-4)
Lifting and holding a frying pan, n=11	1 (1-2)	3 (3-4)
Picking a coin from a table, n=8	2 (1.3-3)	3.5 (2.3-4.0)

Table 1: Percieved limitation

Conclusion: Based on preliminary analysis the glove appears to increase hand function and might improve physical function in persons with IBM who experience impaired hand function. However, persons with severe limitation in combination with reduced function in the arm might not benefit of the glove in daily activites. A prospective intervention study is planned to further investigate the usefulness of the glove.

Malin Regardt, Helene Alexanderson, Stephanie Hunn, Lindsay N Alfano, Roland Mischke, Ingrid de Groot, Anneli Dihkan, Therésia Danielsson, Annika Rydgård, Humza A Chaudhry, Lesley-Ann Saketkoo. A poster was presented at The 5th Global Conference on Myositis (GCOM), Pittsburgh 13-16th of March 2024 by Malin Regardt OT, PhD. at Karolinska Institutet. Results to be published.

Summary of evidence for Carbonhand as an assistive device

Several studies have been performed with Carbonhand, including more than 360 study participants with reduced grip strength and/or impaired hand function due to various reasons and diagnoses.

For example; people with Stroke, Spinal Cord Injury, Multiple Sclerosis, Sarcopenia, and Myopathy benefited from using Carbonhand within these trials.

The results from these trials show:

- Manipulation of several objects encountered in ADL was experienced as improved during the use of Carbonhand
- A sustained and strong grip was perceived as the most beneficial effect while using Carbonhand
- The participants experienced an improved grip strength during the use of Carbonhand
- The usability of Carbonhand was rated as high
- Carbonhand is safe to use.

Reference list:

1. Radder B, Prange-Lasonder G, Kottink A, Melendez-Calderon A, Buurke J, Rietman J. Feasibility of a wearable soft-robotic glove to support impaired hand function in stroke patients. *J Rehabil Med* 2018; 50: 598–606.
2. Radder B, Prange-Lasonder G, Kottink A, Holmberg J, Sletta K, Van Dijk M, Meyer T et al. (2018): The effect of a wearable soft-robotic glove on motor function and functional performance of older adults, *Assistive Technology*, DOI: 10.1080/10400435.2018.1453888
3. Osuagwu B.A. et al. (2019) Clinical Trial of the Soft Extra Muscle Glove to Assess Orthotic and Long-Term Functional Gain Following Chronic Incomplete Tetraplegia: Preliminary Functional Results. In: Masia L., Micera S., Akay M., Pons J. (eds) *Converging Clinical and Engineering Research on Neurorehabilitation III*. ICNR 2018. Pp 385–389. doi.org/10.1007/978-3-030-01845-0_77
4. Nilsson M, Fryxell W. A, Wadell C, von Holst H, Borg J. A Helping Hand – On Innovations for Rehabilitation and Assistive Technology. Grip strengthening glove to improve hand function in patients with neuromuscular disorders. Doctoral Thesis. Stockholm: Royal Institute of Technology; 2013. ISBN 978-91-7501-703-7
5. Hermansson L, Gunnarsson M, Lutterman M. Effect of a grip-strengthening glove on persons with reduced hand function. (2014).
6. Palmcrantz S, Plantin J, Borg J. Factors affecting the usability of an assistive soft robotic glove after stroke or multiple sclerosis. Karolinska Institutet, Department of Clinical Sciences, Danderyd Hospital, Division of Rehabilitation Medicine, Stockholm, Sweden. *Journal of Rehabilitation Medicine* (2020). <https://doi.org/10.2340/16501977-2650>
7. Hashida R, Matsuse H, Totsugi M, Omoto M, Shinozaki N, Tanaka J, Nago T. (2016). Study on the Effects of Motor-Assisted Gloves (SEM™ Glove) on Functional Disorders of the hand. *The Kurume Medical Journal* (2018). DOI:10.2739/kurumemedj.MS652007
8. Radder B, Prange-Lasonder GB, Kottink AIR, Holmberg J, Sletta K, Van Dijk M, et al. The effect of a wearable soft-robotic glove on motor function and functional performance of older adults. *Assist Technol*. 2020;32(1):9–15