

Empowering people with
impaired hand function



A Carbonhand® Hero



Long



Medium



Short



Our mission

We empower people with impaired hand function through innovative and lifechanging solutions.

We Increase independence, improve quality of life, and enable a new Standard of Care.

The Carbonhand[®] system

- A 3-finger glove, with state-of-the-art technology using embedded pressure sensors to enhance grip
- Artificial tendons control the force of each finger to improve grip response and strength
- The Glove can be worn to suit different users and situations
- An App allows for programming of multiple user profiles to fit all user needs and situations
- Battery capacity to facilitate daily activities for a whole day
- Carbonhand provides increased independence and quality of life for the user!



Who can benefit from using Carbonhand?

Indication for Use: People with reduced grip strength and/or impaired hand function.

Examples of patient groups:

Stroke	Stroke
Neurological Injuries	Spinal Cord Injury
	Traumatic Brain Injury
	Brachial Plexus Injury
	Peripheral Nerve Injury
Neurological Illness	Amyotrophic Lateral Sclerosis
	Multiple Sclerosis
	Neuromuscular disorders
Age Related Weakness	Hand Osteoarthritis
	Age related Sarcopenia
	Rheumatoid Arthritis
Musculoskeletal disorders	Carpal Tunnel Syndrome
	Epicondylitis
	Other: Myositis



General Guidelines for Identifying Patients

Inclusion criteria

- Individuals suffering from weakness in hand, wrist and forearm
- Individuals with impaired hand function

Exclusion criteria

- Individuals presenting open wounds or swollen hand before disappearance
- Individuals with severe spasticity in the hand, wrist and forearm

Consider

- Range of Motion (active and passive)
- Voluntary movement and functional ability
- Sensory function
- Pain

User examples



Peter

- Spinal Cord Injury
- Limited functionality in right hand
- Carbonhand provides a firm grip and helps him during daily activities



Helen

- Multiple Sclerosis (MS)
- Weakness and a limited possibility to perform ADL
- Carbonhand gives her better grip and endurance during activities of daily living



Scott

- Post-traumatic hand injury
- Reduced muscle function in the right hand
- Carbonhand makes it possible for Scott to continue his active lifestyle.



Mats

- Peripheral nerve injury
- Impaired functionality of right hand
- Carbonhand gives him a chance to continue to work and play with his kids

Improved grip strength, hand function, and quality of life!

Real Life Examples

Bryan

Diagnosis:
Spinal Cord
Injury

Challenge:
Has a weak grip
in his right hand
and can't hold
objects



John

Diagnosis:
Inclusion Body
Myositis

Challenge:
Has a weak
grip and can't
hold objects



Real Life Examples

Daniel

Diagnosis:
ALS

Challenge:
Weak grip and
tremor



Ingrid

Diagnosis:
Multiple
Sclerosis

Challenge:
Weak grip



Two Wearing Options For Personalized Support

- You can choose to select support for either the index and middle fingers, or the middle and ring fingers
- Enables individuals with varying hand impairments to benefit from the glove
- Wearing options allow flexibility to support different activities



Supports Varying Functional Grip Patterns



Carbonhand - how the system works

- The sensors in the glove provide input when pressure is applied
- The power from the motors is transferred through artificial tendons in the glove
- The behavior of the glove is controlled by a sophisticated control system
- The behavior can be customized for each individual user



Carbonhand® – The Glove



1: Pressure sensors

Dynamic pressure sensors that detects when the user grips an object

2: Artificial tendons

Transfers the forces from the power unit to the fingers in a natural grip movement

3: Control pad

Let the user easily select the preferred profile or initiate lock assist, if this is required by the user

Carbonhand® – The system

1. Glove

Soft, lightweight, and flexible glove with pressure sensors, artificial tendons, and control pad

2. Power unit

Includes battery (8 h), power actuators, control system and glove connector

3. Cord

Contains electrical and bowden cables to transfer sensor signals and forces between the glove and power unit

4. Arm straps

Keeps the cord close to the user in order to avoid entanglement



Ease of use



Open design

And made of flexible material
to fit different anatomies and
impaired hand function



Magnetic buckle

Makes it easy to don and
doff the glove



Control pad

Allows selection of user profiles,
initiate lock assist, or calibration

Adaptable Anywhere, Ready Anytime



Three User Profiles

Each profile offers different settings based on the user's support needs and can be easily adjusted on the control pad.

Further customization and additional profiles are available through the Bioservo App.



Precision Task



Dynamic



Strong & Steady

Bioservo App – Personalized setting of functionality

Carbonhand comes with three pre-defined profiles, that can be personally adjusted to fit the user's needs:



Max force

The maximum force provided by Carbonhand



Stickiness

The force required by the user to release the grip




Responsiveness

The speed between sensor input and force output



Sensor synergies


Linking sensors to multiple fingers





Max Force

- Decides how much force is applied by Carbonhand.
- Can be adjusted for each finger individually to compensate for impairments in specific fingers, or combined if needed equally for all fingers.
- Create specific profiles tailored to your activities and timing, for example, less Max Force in the morning and more in the evening.





Stickiness

- The stickiness helps the user to hold on to objects, but requires more effort to release grip
- At the maximum setting, Carbonhand provides maximum holding stability
- Consider low stickiness setting in patients with weak finger extension
- The grip can also be completely locked by activating the Lock Assist function from the Control Pad





Responsiveness


- Users can adjust how quickly the system responds once a sensor is activated
- An increased responsiveness provides a “magnetic” feeling
- Assists individuals with decreased grip strength to initiate a grip






Linking fingers

- Sensors on one finger can be linked to another finger in order to activate also that finger
- The palm sensor can be used to activate any finger or fingers
- This is useful when the user has difficulties initiating a grip



My Activity - Reporting in the APP

- Motivational monitoring and tracking
- Analyze trends and Track activity
- View the information for one day, one week, three months or one year
- Swipe back and forth in time to see how patterns of use differ.




Remote Connectivity

- Allows clinician or Carbonhand representative to securely log into a patient's Bioservo APP remotely
- Allows remote configuration of glove
 - Add or adjust grip profiles
 - Troubleshoot
- Enables clinicians to access grip reporting



The Bioservo App – Web Demonstration



<https://app.bioservo.cloud/>

Several Carrying Options



Carbonhand can be carried in multiple ways, with a foldable design for carrying in a bag, purse, or mounted on a wheelchair

Evidence on Carbonhand as an Assistive Device

In total 367 patients were included in the trials showing the effect of using Carbonhand. People with reduced grip strength and/or impaired hand function due to various reasons and diagnoses were included as study participants. For example; people with Stroke, Spinal Cord Injury, Multiple Sclerosis, Sarcopenia, and Myopathy benefited from using Carbonhand within these trials.

The results show that:

- Manipulation of several objects encountered in ADL was experienced as improved during the use of Carbonhand
- A sustained and strong grip was perceived as the most beneficial effect while using Carbonhand
- The participants experienced an improved grip strength during the use of Carbonhand
- The usability of Carbonhand was rated as high
- Carbonhand is safe to use.

The medical claim based on these studies:

Carbonhand enhances the user's ability to perform activities in daily living.



Usability Trial – Carbonhand an Assistive Device

Presented at Global Conference on Myositis in Pittsburgh, USA
March 2024

Table 2. Measures on difficulty in chosen activities without and with the glove, (0 = "unable to do", 4 = "without any difficulty").

Chosen activity	Without the glove	With the glove
Lift free weights, n=16	2 (1-2)	4 (3.25-4)
Open previously open jars, n=16	1 (1-2)	4 (3-4)
Lift a heavy bag from the floor, n=15	2 (2-2)	4 (4-4)
Hold a frying pan, n=11	1 (1-2)	3 (3-4)
Pick a coin from a table, n=8	2 (1.3-3)	3.5 (2.3-4.0)

Improved physical function using a power enhancing glove in persons with IBM

Malin Regardt^{1,2}, Helene Alexanderson^{1,2}, Stephanie Hunn¹, Lindsay N Altano³, Roland Mischke⁴, Ingrid de Groot⁵, Anneli Dijkman⁶, Therésa Danielsson⁷, Annika Rydgård⁸, Humza A Chaudhry⁹, Lesley-Ann Sakelko¹⁰ On behalf of MIRA Exercise and Rehabilitation Scientific Working Group.

¹Malmö Red Cross, ²Department of Neurobiology, Care Sciences and Society, Division of Occupational Therapy, Karolinska Institutet, Stockholm, Sweden, ³Warner's Health and Allied Health Professionals Theme, Karolinska University Hospital, Stockholm, Sweden, ⁴Department of Medicine, Solna, Division of Rheumatology, Karolinska Institutet, Stockholm, Sweden and Women's Health, ⁵Department of Neurology, NeuroMuscular Division, Washington University School of Medicine, St. Louis, United States, ⁶The Royal Warne Research Institute at Radcliffe Children's Hospital, Centre for Gene Therapy, Columbia, UK, United States, ⁷The Ohio State University College of Medicine, Department of Pediatrics, Columbus, United States, ⁸Patent Research Partner from Germany, The Netherlands, Sweden, ⁹BioServo Technology AB, Stockholm, Sweden, ¹⁰Louisiana State University and Tulane University Schools of Medicine, New Orleans School of Medicine and Research Center, University Medical Center Comprehensive Pulmonary Hypertension Center and Interstitial Lung Disease Clinic, Program, New Orleans, USA.

Conclusion

Based on preliminary analysis the glove appears to increase hand function and might improve physical function in persons with IBM who experience impaired hand function. However, persons with severe limitation in combination with reduced function in the arm might not benefit of the glove in daily activities. A prospective intervention study is planned to further investigate the usefulness of the glove.

Introduction

Individuals with IBM have reduced hand function leading to limitation in daily activities. There is an urgent need to develop therapies and assistive devices to improve every-day function and quality of life for individuals with IBM.

Aim

To investigate if a power enhancing glove is feasible to use in persons with IBM.

Methods

Data were collected during The Myositis Association's (TMA) patient conference in 2023.

Participants: Individuals with IBM self-reported hand weakness.

The glove's assistive force in grip and pincer strength is triggered by an 'intention-detection' logic that reacts to and supports the follow-through of a hand movement initiation by the user (Figure 1).

Measures: Physical function: IBM-Patient Reported Outcome, Upper Extremity Function Scale (IBM-PRO) and IBM-Functional Rating Scale (IBM-FRS)

Health related quality of life: EQ5D VAS

Grip strength: Jamar

Procedure: Participants selected three activities experienced as difficult and important. They performed the selected activities rating their perceived limitation on a 5-point scale (0 = "unable to do", 4 = "without any difficulty") first without the glove and then upon using the glove. The glove was fitted individually to the left or right hand.

After testing the glove, the participants answered a series of open-ended questions regarding their perception of the glove.

Statistics: Wilcoxon signed rank test were performed to analyze differences between perceived limitation with and without the glove using.

Results

The study included 40 individuals with IBM, (Table 1).


The three most selected activities were *lifting free weights, open previously opened jars, lift a heavy bag from the floor, holding a frying pan and picking a coin from a table.*

All activities were perceived easier to perform with the glove (p<0.039) with ratings on perceived limitation reflecting "with little difficulty" (3) to "without difficulty" (4).

Table 1. Demographic data and status on the participants.

Demographic and status	IBM, n=40
Age md (IQR), years	69 (64-74)
Gender male/female (%)	52.3/47.5
IBM FRS (0-40) md (IQR)	20 (14-25)
PRO-UEFS for IBM (0-48) md (IQR)	23 (15.5-33)
EQ5D VAS (0-100) md (IQR)	74.5 (66.25-83.75)
Jamar right hand md (IQR) kilogram	4.3 (1.6-10.7)
Jamar left hand md (IQR) kilogram	3.2 (1.3-9.6)

IBM=Inclusion body myositis, FRS=Functional Rating Scale, md= median, IQR=interquartile range, EQ5D=EuroQol, VAS=Visual analogue scale.




References for Carbonhand as an Assistive Device

1. Radder B, Prange-Lasonder G, Kottink A, Melendez-Calderon A, Buurke J, Rietman J. Feasibility of a wearable soft-robotic glove to support impaired hand function in stroke patients. *J Rehabil Med* 2018; 50: 598–606.
2. Radder B, Prange-Lasonder G, Kottink A, Holmberg J, Sletta K, Van Dijk M, Meyer T et al. (2018): The effect of a wearable soft-robotic glove on motor function and functional performance of older adults, *Assistive Technology*, DOI: 10.1080/10400435.2018.1453888
3. Osuagwu B.A. et al. (2019) Clinical Trial of the Soft Extra Muscle Glove to Assess Orthotic and Long-Term Functional Gain Following Chronic Incomplete Tetraplegia: Preliminary Functional Results. In: Masia L., Micera S., Akay M., Pons J. (eds) *Converging Clinical and Engineering Research on Neurorehabilitation III*. ICNR 2018. Pp 385-389. doi.org/10.1007/978-3-030-01845-0_77
4. Nilsson M, Fryxell W. A, Wadell C, von Holst H, Borg J. A Helping Hand – On Innovations for Rehabilitation and Assistive Technology. Grip strengthening glove to improve hand function in patients with neuromuscular disorders. Doctoral Thesis. Stockholm: Royal Institute of Technology; 2013. ISBN 978-91-7501-703-7
5. Hermansson L, Gunnarsson M, Lutterman M. Effect of a grip-strengthening glove on persons with reduced hand function. (2014).
6. Palmcrantz S, Plantin J, Borg J. Factors affecting the usability of an assistive soft robotic glove after stroke or multiple sclerosis. Karolinska Institutet, Department of Clinical Sciences, Danderyd Hospital, Division of Rehabilitation Medicine, Stockholm, Sweden. *Journal of Rehabilitation Medicine* (2020). <https://doi.org/10.2340/16501977-2650>
7. Hashida R, Matsuse H, Totsugi M, Omoto M, Shinozaki N, Tanaka J, Nago T. (2016). Study on the Effects of Motor-Assisted Gloves (SEM™ Glove) on Functional Disorders of the hand. *The Kurume Medical Journal* (2018). DOI:10.2739/kurumemedj.MS652007
8. Radder B, Prange-Lasonder GB, Kottink AIR, Holmberg J, Sletta K, Van Dijk M, et al. The effect of a wearable soft-robotic glove on motor function and functional performance of older adults. *Assist Technol*. 2020;32(1):9–15

Clinicians: Education and Support

- Carbonhand Academy with Step-by-step video tutorials
- Exam for Carbonhand Certification
- A Center of Excellence program is available including:
 - Product Training
 - Clinical Training
 - Peer-to-peer Consultations
 - Support



The screenshot shows the Bioservo Academy Carbonhand website. The header includes the Bioservo logo and navigation links: Products, Corporate, News & Media, Bioservo Academy, and Contact. A large banner image features a woman sitting on a couch, holding a cup. The main heading reads "Bioservo Academy Carbonhand". Below the banner, a welcome message states: "Welcome to Bioservo Academy, the portal where you will find all educational training material and documentation for the Carbonhand product. All in one place." A table lists various sections with progress bars and thumbs-up icons. To the right, there is a "Sign out" link and a language selection menu for Carbonhand 2.0, offering options in English, Deutsch, Français, Italiano, Norsk, and Svenska. At the bottom, a call to action encourages users to "Take the exam and get your personalized certificate as proof of your Carbonhand knowledge" and provides a "Begin the examination" button.

Section	New content	Progress
Introduction to the Carbonhand Academy		👍
Introducing Carbonhand		👍
Wearing and handling the system		👍
Glove configuration		
Examples of activities of daily living		
Clinical evidence of using Carbonhand		
Support and documentation		👍

Take the exam and get your personalized certificate as proof of your Carbonhand knowledge
✍️ Begin the examination



Carbonhand – The Components 1/2

Glove with Glove Connector



Power Unit



Shoulder Straps



Arm straps (3)
2-S/M 1-L



Power Unit Sleeve
(Single or Double Carry)



Belt
(S/L, Med=S+L)



Carbonhand – The Components 2/2

Charger



Soft Case



Padded Washing bag



External Button
(optional)



Wheelchair Attachment Straps
(optional)



Power transmission



Ball screws

Transfers the rotation from the electrical motors to a linear movement.

Artificial tendons

in the glove connection are pulled with a maximum force of 20 N each.